


## STEFAN'S CONSTANT MODEL: VTPL-SC-20



- To determine the value of Stefan's constant by measuring temperature difference through thermocouple. The experimental setup should include
- Stefan's apparatus:- black body and thermocouples
- Thermometers:- mercury, 0-360 degree celcius, 3 no.'s
- Stop watch
- Sensitive galvanometer
- Beaker:- Glass
- Connecting Wires
- Steam Boiler
- Electrical heating elements (Hot plates)

*Note: There may be any change in specification due to continuous R & D without notice.*

**VIJAYANTA TECHNOLOGIES PVT. LTD.**

*(Formerly Vijai Electronics)*

Dr. Baldev Singh Marg 28/147 Civil Lines, Roorkee-247667 Distt. Haridwar, Uttarakhand

Phone No.: 01332 – 272509, 7579200827

E-Mail : vijayantatechologies@gmail.com, vijaielectronics1965@gmail.com